

# Product Catalogue



Mission  
To Care  
To Cure



## HOSPITAL CLOTHING

### Linen Medical Clothing

- Surgeon and medical staff clothing/Apparel
- Patient clothing
- Hospital bedding and draping

### Disposable Medical Range

- Disposable Surgeon and medical staff clothing/Apparel
- Disposable Patient clothing
- Disposable Hospital bedding, draping & medical kit



We are manufacturer and supplier of an extensive range of linen medical cloth used in hospitals and various medical centers. Our all products are made in the supervision of experience professionals which works on only one motto  
“SAFETY AND PROTECTION OF OUR CUSTOMER FIRST”

The fabric used in our products are eco-friendly and reliable to use.



### Laboratory Coat

We offer a superior range of fabric lab coat to our clients at competitive price.

Size : S, M, L, XL, XXL

Material used : Cotton (tusser and spun)

Color : White/Sky blue/Green/Grey & Navy blue.

**(Available in half/Full sleeves with three pockets)**

### Surgeon's Gown

We offer a superior quality of surgeons gown specially designed for doctors and surgeons to be worn at the time of surgical procedures for their safety and protection and prevent them from medium level splashes.

Size : Standard

Material used : Cotton casement

Color : Green and sky blue



### Surgeon's Suit

We offer a superior quality surgeons suit which provides them comfort dryness and activeness

Size : Standard

Material used : Cotton

Color : Pink/Sky blue/Green/Maroon

## Patient Suit

We offer a superior quality patient suit which provides them comfort dryness and activeness

**Size :** Standard

**Material used :** Micro Light Weight Cotton

**Color :** Sky blue/Green/Mehndi



## Patient Gown

We are engaged in offering a qualitative range of patient gown

**Size :** Standard

**Material used :** Cotton century

**Color :** White/blue check

## Patient Gown with Pajama

We offer a wide range of patient gown with pajama for the comfort of our clients

**Size :** Standard

**Material used :** Cotton

**Color :** White/Grey lining



## Surgeon's Cap & Mask

With rich industry experience, we are manufacturing and supplying superior surgeons cap and mask.

**Size :** Standard

**Material used :** Casement cotton

**Color :** Green/Blue

## Plastic Apron (auto cleavable)

We offer a superior quality plastic apron which are auto cleavable and supply after undergo various quality parameters

**Size :** Standard

**Color :** Navy blue





## Hospital Bed Sheets

We provide a wide range of Bed sheets to our clients at industries best price.

Size : 150 cm X 250 cm

Material used : Casement heavy quality and dasooti

Color : Violet/Indigo/Blue/Green/Yellow/Orange/Red/White/Sky blue


## Hospital Towels

We provide a wide range of cotton Towels to our clients at industries best price.

Size : Depend on order

Material used : Cotton

Color : White/Green/Sky blue/Dark green



## Hospital Blankets

We are providing multi-color hospital blankets of superior quality to our clients at industries best price.

Size : 225 X 150 cm

Material used : Cotton/Rashmi

Weight : 2 kg to 2.5 kg

Color : Dark green/Baze/Orange/Red/Navy blue

## O.T./ Draw Sheets (Fenestrated & Non-fenestrated)

We provide a wide range of O.T./ draw sheets to our customer at industries best price. These sheets are available in fenestrated as well as non-fenestrated.

Size : 90 cm X 90 cm, 90 cm X 145 cm, 145 cm X 180 cm, 150 cm X 220 cm

Material used : Casement heavy quality

Color : Green/White/Sky blue



## Mackintosh (Auto-cleavable)

We provide a wide range of waterproof Mackintosh to our clients at industries best price.

Size : 100 X 150 cm,

Material used : Matte

Color : Black/Navy blue



### Disposable Surgeon's Gown (Reinforced)

We are engaged in offering high quality medical grade non-woven Disposable Surgeon's Gown which is separately sterilized by ETO gas to provide 100 % bacteria free environment surgeon.

Size : Standard

Material Used : SMMMS non-woven fabric (Anti-static)

Colour : Blue

(Reinforced gown is having special layer on front-side and inside arms)

### Disposable Surgeon's Gown (Regular)

We are engaged in offering high quality medical grade non-woven Disposable Surgeon's Gown which is separately sterilized by ETO gas to provide 100 % bacteria free environment surgeon.

Size : Standard

Material Used : SMMMS non-woven fabric (Anti-static)

Colour : Blue



### Disposable Wrap Round Surgeon's Gown (Full Gown)

We are engaged in offering high quality medical grade non-woven Disposable Wrap Round Surgeon's Gown which is separately sterilized by ETO gas to provide 100 % bacteria free environment surgeon.

Size : Standard

Material Used : SMMMS non-woven fabric (Anti-static)

Colour : Blue

### Disposable Face Mask/ Cramp cap/ Surgeon cap/ shoe cover

Made from high quality medical grade non-woven cloth


- 3 ply mask tie type/elastic loop in single packing sterilized
- Surgeon Cap in single packing sterilized
- 2 ply/3 ply Mask/surgeon cap/cramp cap

(Available in bulk packing too)



### Disposable Patient suit

Made from high quality medical grade non-woven cloth  
Each patient suit pack sterilize by ETO gas.



## HIV/AIDS/Hepatitis prevention Kit

This pack contains :

- ⇒ Disposable Surgeon gown (SMMMS-Reinforced Gown)
- ⇒ Disposable Face Mask (3 Ply with bacteria filter)
- ⇒ Disposable Cap
- ⇒ Disposable Bed Sheet
- ⇒ Disposable Shoe Cover (Long)
- ⇒ Surgical Gloves (Elbow Length)
- ⇒ Goggles
- ⇒ Waste Carry Bag
- ⇒ Towel

(This pack individually sterilized by ETO gas)

## Disposable Delivery Kit

This pack contains :

- ⇒ Disposable surgeon gown (regular)
- ⇒ Dripping Sheet
- ⇒ Baby wrap towel (60 cm X 80 cm)
- ⇒ Surgical gloves (2 pairs)
- ⇒ Disposable face mask
- ⇒ Disposable shoe cover
- ⇒ Surgical blade (11 no.)
- ⇒ Disposable Bed sheet PVC (120 cm X 150 cm)
- ⇒ Disposable PVC apron, Cord clamp
- ⇒ Gauze swab (4 PCs), Cramp cap (2 Pcs), Paper Napkin

(This pack individually sterilized by ETO gas)



## General Surgery Kit/ O.T. Kit

This pack contains :

- ⇒ Disposable surgeon gown (reinforced)
- ⇒ Disposable PVC bed sheet (120 cm X 210 cm)
- ⇒ Disposable surgeon gown
- ⇒ Surgical glove (No. 7.5)
- ⇒ Surgical glove (No. 8)
- ⇒ Disposable face mask
- ⇒ Disposable shoe cover (Non-woven fabric)
- ⇒ Disposable towel

(This pack individually sterilized by ETO gas)

## Other Disposable Range



Disposable Draw Sheet  
(Non-woven fabric)



Disposable Bed Sheet  
(PVC/PE Material)



Baby Wrap Towel  
(White Spun Bond)



Half Gown  
(PE Polythene)



Dental Apron  
(Laminated White Spun Bond)



Disposable O.T. Sponge (Mop)  
(Cotton with X-ray Thread)

**Exclusive Range : C-arm Cover, Fluoroscope Cover, PCNL Drape, Angiography Drape**



A Quality Product from :  
**Perfect Surgicare Industries Pvt. Ltd.**  
 Plot no.-11, Ground Floor, Vijay Nagar,  
 Near Namastey Chowk, KARNAL-132001  
 (HARYANA) INDIA  
 E-mail : Perfectsurgicare123@gmail.com  
 Customer Care No. : +91-9416100857  
 Mfd. by : Perfect Surgical Industries, Karnal (Hr.)

ISO 9001 : 2015 G.M.P. CE  
 CERTIFIED COMPANY

